

CD45RO (T-Cell Marker) Ab-1

Catalog # MS-112-P0, -P1, or -P (0.1ml, 0.5ml, or 1.0ml at 200ug/ml)

Catalog # MS-112-R7 (7.0ml)

Catalog # MS-112-PCS

INTENDED USE:

- **For In Vitro Diagnostic Use:** This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy.
- **Description:** CD45RO (an isoform of leukocyte common antigen) is expressed on most thymocytes, monocytes, macrophages and granulocytes. It is not expressed on normal B lymphocytes. Intrathymic T cell maturation progresses from CD45RA-/CD45RO+ to CD45RA+/CD45RO-; while peripheral T cell maturation progresses from CD45RA+/CD45RO- to CD45RA-/CD45RO+.
- **Expected Staining Pattern:** Cell Membrane
- **Positive Control:** Tonsil or lymph node

MATERIALS PROVIDED:

CD45RO (T-Cell Marker) Ab-1 (refer to catalog number):

- #MS-112-P (or -P0, -P1): 200ug/ml of antibody purified from ascites. Prepared in 10mM PBS, pH 7.4, with 0.2% BSA and 0.09% sodium azide.
or
- #MS-112-R7: (7.0ml) of antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing stabilizing protein and 0.015mol/L sodium azide.
or
- #MS-112-PCS: 5 positive control slides.
- **Antibody Concentration:** 200ug/ml
- **Host:** Mouse
- **Mol. Wt. of Antigen:** 180kDa
- **Epitope:** Not determined
- **Species Reactivity:** Human, Rhesus Monkey, and Chimpanzee. Does not react with rat. Others-not known.
- **Clone Designation:** UCHL-1
- **Ig Isotype / Light Chain:** IgG2a / kappa
- **Immunogen:** Interleukin-2-dependent human T lymphocytes
- **Microbiological State:** This product is not sterile.

MATERIALS REQUIRED, BUT NOT PROVIDED:

- **Antibody Diluent:** For concentrated antibodies, the antibody must be diluted before using. Use Lab Vision Antibody Diluent (catalog # TA-125-UD). Refer to diluent product instructions for use.
- **Negative Control Reagent:** Refer to the "General Protocol" instructions.
- **Visualization System:** Refer to the "General Protocol" instructions.

METHODS AND PROCEDURES:

Specimen Preparation	Refer to the "General Protocol" instructions.
Dilution of Concentrated Antibody	1:100-1:200 in antibody diluent
Tissue Section Pretreatment	No special pretreatment is required for immuno-histochemical staining of formalin-fixed tissues.
Primary Antibody Incubation Time	30 minutes at Room Temperature
Visualization	To detect antibody, follow the instructions provided with the visualization system.

STORAGE and STABILITY:

This product contains sodium azide and is stable for 24 months when stored at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent is not stored as recommended, performance must be validated by the user.

CD45RO (T-Cell Marker) Ab-1

Catalog # MS-112-P0, -P1, or -P (0.1ml, 0.5ml, or 1.0ml at 200ug/ml)

Catalog # MS-112-R7 (7.0ml)

Catalog # MS-112-PCS


REFERENCES:

- 1) Smith S H, et. al. (1986) Immunology, 58:63-70.
- 2) Leukocyte Typing III, p225-227.
- 3) Akbar A. N, et. al. (1988) J Immunol, 140:2171-2178.
- 4) Nature, 360:264-265.

Thermo Fisher Scientific
Anatomical Pathology
47777 Warm Springs Blvd.
Fremont, CA 94539, USA
Tel: 1-510-991-2800
Fax: 1-510-991-2826
<http://www.thermo.com/labvision>



Manufactured by:
NeoMarkers
For
Lab Vision Corporation



Thermo Fisher Scientific
Anatomical Pathology
93-96 Chadwick Road, Astmoor
Runcorn, Cheshire WA7 1PR, UK
Tel: 44-1928-562600
Fax: 44-1928-562627
Labvision.uk@thermofisher.com