

## WT1 clone WT49


### Instructions for Use

#### Specification:

The WT1 gene, identified as a tumor suppressor gene located at 11p13, is involved in the development of Wilms' tumor. The WT1 gene encodes a transcription factor with four DNA-binding zinc fingers at the C terminus. In vitro studies showed that WT1 suppresses or activates a number of genes, including those for PDGF-A chain, EGF receptor, CSF-1, IGF-II, IGF-I receptor, RAR- $\alpha$ , c-myc, bcl-2 and WT1 itself. Immunohistochemically, WT1 is detected in the nucleus of tumor cells of Wilms' tumor and mesothelioma; therefore, WT1 has traditionally been used as a diagnostic marker for these tumors. Recent reports showed that other types of cancers, such as ovarian serous cancers and rhabdomyosarcomas also express WT1.

#### Availability:

Catalog No.	Contents	Volume
ILM3951-C01	WT1 clone WT49	0,1 ml concentrate
ILM3951-C05	WT1 clone WT49	0,5 ml concentrate
ILM3951-C1	WT1 clone WT49	1,0 ml concentrate



**Reactivity:** Human

**Clone:** WT49

**Species of origin:** Mouse

**Isotype:** IgG<sub>1</sub>/k

**Control Tissue:** Ovarian carcinoma (non-mucinous carcinoma), malignant mesothelioma, testes, kidney

**Staining:** Nuclear

**Presentation:** Tissue culture supernatant containing 15mM sodium azide

#### Application and suggested dilutions:

Pre-treatment: Heat induced epitope retrieval in 50 mM Tris buffer pH9.5, for 20 minutes is required for IHC staining on formalin-fixed, paraffin embedded tissue sections.

- Immunohistochemical staining of formalin-fixed, paraffin embedded tissue section (dilution up to 1:800 / 1:1600)

The optimal dilution for a specific application should be determined by the investigator.

**Note:** Dilution of the antibody in 10% normal goat serum followed by a goat anti-mouse secondary antibody-based detection is recommended.

**Storage & Stability:** Store at 2-8 °C. Do not use after expiration date printed on the vial.

#### References:

- 1) Nakatsuka A et al, Modern pathology, 2006, 19(6): 804-814
- 2) Acs G et al, International Journal of Gynecologic Pathology, 2004, 23(2): 110-118
- 3) Goldstein N S et al, American journal of clinical pathology, 2002, 117(4): 541-545
- 4) Grubb G R et al, Laboratory investigation; a journal of technical methods and pathology, 1994, 71(4): 472-479
- 5) Coosemans A, et al, Gynecologic oncology, 2008, 111(3): 502-508