

UltraMount Plus

INTENDED USE

For In Vitro Diagnostic Use

<u>AVAILABILITY:</u>	<u>Catalog #</u>	<u>Volume</u>
	TA-060-UMX	60 ml
	TA-125-UMX	125 ml

<u>SPECIFICITY:</u>	N/A
<u>ENZYME:</u>	N/A
<u>CHROMOGEN/SUBSTRATE:</u>	N/A

DESCRIPTION

UltraMount Plus is formulated for covering tissue sections and cell preparations immunohistochemically stained with chromogens producing alcohol soluble OR alcohol insoluble end products.

UltraMount Plus can be used with chromogens such as Fast Red or AEC that require an aqueous mounting medium, or those that do not, such as DAB. No heating is required prior to use. This product is optimal for use without cover slips.

WARNINGS & PRECAUTIONS

Refer to MSDS.

STORAGE & SHELF LIFE

Store at room temperature. Do not use past the expiration date.

MICROBIOLOGICAL STATE

Product(s) not sterile.

MATERIALS REQUIRED BUT NOT PROVIDED

N/A

SPECIMEN & REAGENT PREPARATION

Refer to Procedure.

PROCEDURE

Apply several drops of UltraMount to stained specimens. Allow to dry.

REFERENCES

N/A

TROUBLESHOOTING

Please contact Thermo Fisher Scientific Technical Support by phone (1-510-991-2800 or 1-800-828-1628) or by email (lab.reagents@thermofisher.com).

